



TECHNICAL TIP

OVERVIEW OF A PARAMETRIC RELEASE VALIDATION FOR EO STERILIZATION

Background: Traditional EO Monitoring

Traditionally, ethylene oxide (EO) sterilization processes are routinely monitored with biological indicators composed of the microorganism *Bacillus atrophaeus*, usually inoculated onto a 1.5 x 0.25 inch paper carrier. A good example of a biological indicator is the Spordex® (STERIS) bacterial test strip.

For routine monitoring purposes, the biological indicators are placed throughout the sterilizer load and subjected to the sterilization process. After the process, the biological indicators are removed from the load and forwarded to the testing laboratory, where they are placed into a special growth medium and subjected to ideal growth conditions for 7 days (as addressed in USP1). After the 7-day incubation period, negative growth of the biological indicators demonstrates that the sterilization process was effective. The sterilizer load may then be considered for release to market, provided all other release criteria are met.

When using this traditional biological indicator monitoring of the EO sterilization process, expenses are incurred in two ways. First, there is the expense of the biological indicators and the associated laboratory testing. Secondly, there is a large amount of capital tied up in the inventory which must be held “on status” until after the laboratory testing is complete (7 days).

In addition to the expense, the use of biological indicators will increase risks due to the possibility of contamination of the indicator while being handled at the testing laboratory. In most cases, a failure of a biological indicator, even though laboratory induced, would be considered a sterility failure and result in the requirement for reprocessing of all of the materials which were contained in the sterilizer load. This increases costs both of the monitoring system and held inventory.

An Alternative Monitoring Method

ANSI/AAMI/ISO 11135:2014, titled “Sterilization of health care products - Ethylene oxide -: Requirements for the development, validation, and routine control of a sterilization process for medical devices” defines Parametric Release as “declaration that product

is sterile, based on records demonstrating that the process parameters were delivered within specified tolerances.”

Simply, parametric release allows product to be released to the market based only on process records instead of the traditional biological indicator sterility test. This is advantageous in that it eliminates the routine costs and risks associated with biological indicators and laboratory testing. In addition, most companies realize additional savings associated with a reduced unreleased inventory time (providing that the EO residue hold time is shorter than the biological incubation times).

Process Performance Qualification (Microbiological)

ANSI/AAMI/ISO 11135 provides options for the process performance qualification of the EO sterilization process as reference Annex A and Annex B of the standard (note: Additional guidance may be found in ISO 14161).

Annex A, titled “Determination of lethal rate of the sterilization process – Biological indicator/bioburden approach” outlines a method where resistance determinations are demonstrated for the biological indicator and compared to the natural product bioburden.

This is accomplished by running a predefined process at graded exposure times and determining the lethal rate (rate of inactivation) delivered. The knowledge of this rate and the population and relative resistance of the bioburden allows one to establish an exposure time so that a Sterility Assurance Level (SAL) can be predicted.

Annex A outlines two options for data analysis.

1. Direct enumeration: Direct enumeration is the process of determining the lethality of the sterilization process by construction of a survivor curve using direct enumeration (physical counts through serial dilutions) of surviving organisms. At least five points employing graded exposure times, with all other parameters (except time) remaining constant, are utilized. The data generated will enable the

FOR MORE INFORMATION

STERIS Applied Sterilization Technologies

Web: www.steris-ast.com // Email: ast_info@steris.com

(EMEA) +44 (0) 8456 88 99 70

(Americas) 877.783.7479





TECHNICAL TIP

OVERVIEW OF A PARAMETRIC RELEASE VALIDATION FOR EO STERILIZATION

calculation of the time of exposure needed to achieve sterility of the biological indicator.

2. Fraction negative: The fraction negative method also requires graded exposure times to assure survivors, but the post-processing testing methods are different. For this method, a minimum of five exposure times are required. After exposure to the process, the samples are assayed by direct immersion into the appropriate culture medium (pass/fail) in lieu of the physical count performed in Method A. Using the data generated and statistical models provided in the standard (Holcomb-Spearman Karber or Stumbo Murphy Cochran), the death kinetics or D-value, may be calculated. Using the D-value data, an exposure time needed to achieve the desired sterility assurance level can be determined.

TechTeam Discussion

Historically, the use of “process development” by direct enumeration or fraction negative methods (D-value determination) were requirements (ANSI/AAMI/ISO 11135:1994) for validation of a parametric release process. It was the opinion of industry experts (ISO) that these validation methods provided a better understanding of the lethality delivered by the chosen process thus should be a requirement for parametric release.

Unfortunately, the requirements for graded exposure times which promote partial survivors of the biological indicators are extremely difficult to achieve and, in many cases, lead to an extensive and expensive process development/validation effort. The complexity of the program and financial burden associated with the expensive laboratory testing prevented many clients from pursuing the parametric release options.

After many discussions within the ISO group, the 2007 issue of the standard was revised to allow validation of a parametric release process using the overkill method. The overkill method has been in use for many years and continues to be the predominant method of EO Process Performance Qualification within the industry today. Annex B identifies the methods utilized to validate with the overkill method.

Annex B: “Conservative determination of lethal rate of the sterilization process – Overkill approach”.

The overkill method requires a total of three consecutive (1/2 EO exposure time) cycles be performed which result in total inactivation of the biological indicators (of which the initial population was not less than 10⁶). By demonstrating the inactivation of the 10⁶ BI using one half of the exposure time, a Sterility Assurance Level (SAL) 10⁻⁶ is assured when the exposure time is doubled for the routine full cycle.

In addition to the three successful half cycles, the standard requires a cycle of short duration (fractional) from which survivors can be recovered be performed to demonstrate (validate) the adequacy of the recovery technique. Also, it is during this fractional cycle that the resistance of the bioburden is proven to be equal to or less than the resistance of the biological indicator.

Process Performance Qualification (Physical)

As outlined in the standard, the physical performance qualification shall demonstrate:

1. Reproducibility of the process, and shall include a minimum of three consecutive, planned qualification runs in which all the specified criteria are met;
2. That the specified acceptance criteria are met throughout the load for the duration of the proposed routine specification.

Elements of the physical PQ may be conducted during the microbiological PQ. If a) is performed in parallel with the microbiological PQ, then at least one additional qualification run shall be performed in compliance with this requirement.

The physical PQ shall confirm the process such that:

1. The minimum temperature of product to enter the sterilization process and/or the defined conditions required to achieve it shall be established;
2. At the end of the defined preconditioning time (if used), the sterilization load temperature and humidity have been established;

FOR MORE INFORMATION

STERIS Applied Sterilization Technologies

Web: www.steris-ast.com // Email: ast_info@steris.com

(EMEA) +44 (0) 8456 88 99 70

(Americas) 877.783.7479





TECHNICAL TIP

OVERVIEW OF A PARAMETRIC RELEASE VALIDATION FOR EO STERILIZATION

3. The specified maximum elapse time between the completion of preconditioning (if used) and the commencement of the sterilization cycle is appropriate;
4. At the end of the defined conditioning time, if used, the sterilization load temperature and humidity have been established;
5. The chamber humidity was recorded if parametric release was to be used;
6. Gaseous ethylene oxide has been admitted to the sterilizer chamber;
7. Pressure rise and the quantity of ethylene oxide used or the concentration of ethylene oxide in the sterilizer chamber have been established;
8. During the sterilization cycle, the temperature and humidity (if recorded) of the chamber and, where applicable, other process parameters have been established;
9. The temperature of the product load during exposure has been established;
10. During aeration, the temperature of the sterilization load has been established.

For establishments that have widely varying load configurations, the extent to which the variation affects the sterilization process shall be evaluated. It shall be demonstrated that all product exposed to a sterilization process shall be evaluated. It shall be demonstrated that all product exposed to a sterilization process achieves the required SAL.

Review and Approval of the Validation

Upon completion of the validation effort, a report shall be prepared which describes or references specific validated product, the defined loading patterns, and the documented specification for the process. The report will include the value and tolerances for:

1. The minimum temperature of product prior to enter the

sterilization process and/or the defined conditions required to achieve the minimum required;

2. Preconditioning (if used):
 - a. Time in chamber/area, temperature and humidity of chamber/area;
 - b. Temperature and humidity of the sterilization load;
 - c. Maximum elapse time between removal of the load from preconditioning time and commencement of the sterilization cycle;
3. Vacuum levels and rate of evacuation (If used):
 - a. Holding time under vacuum (if used);
4. Inert gas flushing:
 - a. Pressure (change or terminal pressure) and rate (P/ time) of attainment of pressure associated with inert gas/ steam;
 - b. Depth (P or terminal pressure) and rate (P/time) of attainment of vacuum;
 - c. Number of times of repetition and any variations in successive repetitions;
5. Conditioning and/or humidity dwell phase (if used):
 - a. Pressure levels and/or rate of attainment of vacuum or relative humidity levels (whichever is being controlled and monitored);
 - b. Number of steam pulses/vacuum (if used);
 - c. Time;
 - d. Chamber temperature;
 - e. Temperature and humidity of the sterilization load at the end of conditioning;
6. EO Injection and exposure:

FOR MORE INFORMATION

STERIS Applied Sterilization Technologies

Web: www.steris-ast.com // Email: ast_info@steris.com

(EMEA) +44 (0) 8456 88 99 70

(Americas) 877.783.7479





TECHNICAL TIP

OVERVIEW OF A PARAMETRIC RELEASE VALIDATION FOR EO STERILIZATION

- a. EO injection pressure rise (P), EO injection time and terminal pressure of EO injection phase;
 - b. Evidence that the gaseous EO has been admitted to the sterilization chamber by the pressure rise and by one of the following:
 - i. Mass of EO used;
 - ii. Volume of EO used;
 - iii. Direct measurement of EO concentration;
 - iv. Volume of EO used.
 - c. Sterilizer chamber temperature;
 - d. Exposure time;
 - e. Temperature of the sterilization load;
 - f. An indication of the satisfactory operation of the chamber gas circulation system (if used) during exposure;
7. Post exposure flushing (if used):
- a. Depth (P or terminal pressure) and rate (P /time) of attainment of vacuum;
 - b. Pressure (P or terminal pressure) and rate (P /time) of attainment of pressure associated with inert gas/air/steam;
 - c. Number of times of repetition and any variations in successive repetitions;
8. Aeration (if used): Time and temperature;
- a. Time and temperature within the chamber and/or room;
 - b. Pressure changes (if any) in the chamber and/or room;
 - c. Rate of change of air or other gas;
 - d. Temperature of the sterilizer load.

If parametric release is to be used, the validation shall establish:

1. The value and tolerances for chamber humidity by direct measurement during conditioning;
2. The value and tolerance for the ethylene oxide concentration, determined from direct analysis of the chamber atmosphere using analytical methods to establish the process specification for routine processing. The sampling shall be conducted at defined intervals sufficient to verify the required conditions throughout EO Exposure.
3. Temperature of the chamber; recorded from two separate monitoring locations.

Requalification Requirements

Per ISO 11135:2014 the following requirements (the same as a biological release validation) shall apply:

1. Requalification shall be reviewed annually to determine the extent of requalification that is necessary. This shall include an assessment of the need to reconfirm the product SAL through microbiological studies. The outcomes of this review, including the rationale for decisions reached, shall be documented;
2. Requalification of a sterilization process carried out with specified equipment shall be performed at defined intervals against specified acceptance criteria and in accordance with documented procedures. These intervals shall be justified.
3. If requalification indicates that the sterilization process might no longer be capable of achieving the required product SAL, the cause shall be investigated and corrective and/or preventive action shall be taken. As part of the investigation, the effect on the achievement of the specified SAL for previously processed loads of product shall be considered and a risk assessment undertaken on their suitability for use. If investigation shows that the required SAL can no longer be achieved, then a new MPQ/PPQ shall be performed to re-establish the required

FOR MORE INFORMATION

STERIS Applied Sterilization Technologies
Web: www.steris-ast.com // Email: ast_info@steris.com
(EMEA) +44 (0) 8456 88 99 70
(Americas) 877.783.7479





TECHNICAL TIP

OVERVIEW OF A PARAMETRIC RELEASE VALIDATION FOR EO STERILIZATION

SAL. The investigation and subsequent actions shall be record.

Typically, a re-qualification cycle (Half Cycle) is run every two years and a documented review is performed in those years when the re-qualification cycle is not executed. However, depending on the results of the documented review, further work might be required.

Special Considerations for Parametric Release

When designing a parametric release program, special adjustments are made to the validation parameters to assure that the “routine process” falls within “acceptable limits” each time a load is routinely processed. This is accomplished by validating at “worst case” conditions designed to provide a wide “processing window” for routine processing.

1. Half cycles: Simulating winter conditions for all half cycles: To assure that preconditioning (if used) is effective regardless of climatic conditions or facility location, each half cycle validation load is exposed to simulated winter conditions prior to half cycle processing. Successful biological results when challenged with this worst case load conditioning assures that the process is effective for all conditions regardless of the time of the year or the location of the facility.

Relative humidity during conditioning: During each half cycle, the set point for conditioning is set below the nominal (routine) value to demonstrate adequate process lethality in the lower region of the relative humidity processing range. The relative humidity value reported for the successful half cycle is used as the routine minimum for the end of conditioning which is a requirement for parametric release.

Establishing a minimum EO concentration during exposure: During each half cycle, the set point for gas injection is reduced compared to the nominal (routine) value, at least 50 mg/L. The minimum concentration reported for the successful half cycle is the routine minimum value for exposure which is a

requirement for parametric release.

2. Full cycles: During each full cycle (full loads), the set point for gas injection is set for the nominal (routine) value plus 50 mg/L or more. The maximum concentration reported for the
3. Mixed loads: For those Customers that present mixed loads to the sterilizing facility, a fourth half cycle is recommended using the minimum loading configuration at minimum density. This load is processed in a biological challenged half cycle using minimum cycle conditions. Successful biological results using this minimum loading condition when supported by the successful maximum loading half cycles, serves to validate the process for any load mix from the minimum loading qualified to the maximum loading qualified.

FOR MORE INFORMATION

STERIS Applied Sterilization Technologies
Web: www.steris-ast.com // Email: ast_info@steris.com
(EMEA) +44 (0) 8456 88 99 70
(Americas) 877.783.7479





TECHNICAL TIP

OVERVIEW OF A PARAMETRIC RELEASE VALIDATION FOR EO STERILIZATION

REFERENCES

1. United States Pharmacopeia, 31st edition, 2008.
2. ANSI/AAMI/ISO 11135:2014, Sterilization of health care products – Ethylene oxide: Requirements for the development, validation, and routine control of a sterilization process for medical devices

FOR MORE INFORMATION

STERIS Applied Sterilization Technologies
Web: www.steris-ast.com // Email: ast_info@steris.com
(EMEA) +44 (0) 8456 88 99 70
(Americas) 877.783.7479

